

Use this document to log each of your QC check records.

### From the MOT Testing Guide (Section B6 - Assurance)

The frequency of checks may typically be expected to be 1 per tester every 2 months, however this is based on an average throughput of only 2-3 tests per day (for experienced testers). It is expected that as the volume of MOT tests increases, so too should the frequency of QC checks. MOT Juice recommends a minimum of 1 QC check per tester per month.

*Please note: This log alone may not be sufficient evidence of QC checks being carried out, and you should aim to retain as many backup documents to prove your QC checks took place.*

Name: \_\_\_\_\_

Date of check	Checked by	QC Comments

### Need more rows?

Download additional pages here: <https://www.motjuice.co.uk/downloads/qc-log-extension.pdf>

### Did you know?

You can log all your QC checks online with the MOT Juice Quality Control system, so you don't need to retain any paper copies. Your records are stored securely online and ready to print at the click of a button for when a VE visits.